Immersion Swaddle Bathing



Transform bathing to a positive experience



Benefits of Immersion Swaddle bathing:

- Minimize temperature loss 1-3
- Decreases crying 1,3
- Reduces physical and behavioral stress 2, 4-6
- Supports family centered care 6,7,10
- Enhances ability to feed after bath 6,10
- Evidence-Based Practice 1-6, 8-10

Swaddle bathing with the TurtleTub complements delayed bathing and minimizes adverse clinical outcomes associated with traditional bathing methods. After participating in hospital bathing, parents can continue to swaddle bathe at home. Preterm infants, full-term infants, and infants with NAS all benefit from immersion swaddle bathing. In addition, the TurtleTub is a strong parent pleaser, providing a marketing advantage for your hospital.



www.catapult-products.com

Early Experiences Matter

The TurtleTub System:

Designed for infant safety, comfort, ease of use and infection control.



TurtleTub Infant Bathtubs:

- Ideal for swaddle bathing and immersion tub bathing
- Options for single-patient and multi-patient use
- · Ribs and high sides provide excellent support
- Integrated temperature indicator
- Smooth contours for easy cleaning
- Convenient pour spot
- Stackable
- Complies with ASTM F2670 standard for infant bathtubs

TurtleTub Dimensions:

22"L x 18" W x 6"H

TurtleTub Disposable Bath Liners (for multi-patient tub use):

- Strong/thick poly with elastic holds liner securely in TurtleTub
- Infection control compliant
- Non-leaching poly material complies to USP Physiochemical 661 and Cytotoxicity tests providing a safe bath
- · Available in packs with accessories or packaged individually





TurtleTub Bathing Cart:

- 4-sided guard rail keeps TurtleTub secure
- · Two locking wheels for safe bathing
- Two smooth gliding "soft-close" drawers

Specifications:

- 20" wide x 16" deep x 34" height (not including guard rail)
- 4-sided guard rail (adds 3.5" to height)
- Stainless steel grade 304







Optional Cart Caddy (sold seperately)

Cart Caddy hangs conviently on the guard rail to hold bathing supplies

TurtleTub Accessories:

Available individually and in our kits and packs (See ordering info for details)

Baby Wash:

- Single-use packets (10mL)
- pH nuetral and mild
- Made without dyes, sulfates, parabens or phthalates
- Fragrance-free

22 oz. rinse cup











Brush and comb set for full-term infants

Mesh storage bag to hang, dry and store the TurtleTub

TURTLETUB™ ORDERING INFORMATION

Item #	Description	Qty/Case
Complete Swaddle Bathing Kits		
TT1100	TurtleTub Standard Kit Each kit contains: 1 TurtleTub with Integrated Temperature Indicator; 1 Neutral pH Mild Baby Wash; 1 Small Fleece Blanket; 1 Rinse Cup; 1 Scalp brush; 6 Washcloths & Instructions for Use.	10 Kits
TT1100MB	TurtleTub Mother Baby Unit Standard Kit Each kit contains: 1 TurtleTub with Integrated Temperature Indicator; 1 Neutral pH Mild Baby Wash; 1 Large Cotton Flannel Blanket; 1 Rinse Cup; 1 Brush & Comb set; 6 Washcloths & Instructions for Use.	10 Kits
TT1106	TurtleTub NICU Starter Kit w/Liner - Each kit contains: 5 TurtleTub Bathtubs (each with Integrated Temperature Indicator & Rinse Cup); 10 Small Fleece Blankets and 25 Accessory Packs w/Liner (TT1107).	1 Kit
TT1113	TurtleTub Deluxe Kit - Each kit contains: 1 TurtleTub with Integrated Temperature Indicator; 1 Neutral pH Mild Baby Wash; 1 Large Fleece Blanket; 1 Rinse Cup; 1 Scalp brush; 1 Mesh Storage Bag; 6 Washcloths & Instructions for Use.	10 Kits
TurtleTub - Bathtubs (Tub only - no accessories)		
TT1104	TurtleTub Bathtub (with Integrated Temperature Indicator)	10
TT1104-2	TurtleTub Bathtub (with Integrated Temperature Indicator)	2
Disposable Bath Liners & Liner Packs		
TT1124	Mother Baby Unit Premium Bath Accessory Pack Each pack contains: 1 Disposable Liner; 1 Disposable Swaddle Bath Blanket; 1 Neutral pH Mild Baby Wash; 1 Brush & Comb set & 6 Washcloths.	50
TT1125	Mother Baby Unit Deluxe Bath Accessory Pack Each pack contains: 1 Disposable Liner; 1 Disposable Swaddle Bath Blanket; 1 Brush & Comb set & 6 Washcloths.	50
TT1122	Mother Baby Unit Standard Bath Accessory Each pack contains: 1 Disposable Liner; 1 Neutral pH Mild Baby Wash; 1 Brush & Comb set & 6 Washcloths.	50
TT1123	Mother Baby Unit Basic Bath Accessory Pack Each pack contains: 1 Disposable Liner; 1 Brush & Comb set & 6 Washcloths.	50
TT1107	NICU Standard Bath Accessory Pack - Each pack contains: 1 Disposable Liner; 1 Neutral pH Mild Baby Wash; 1 Scalp brush & 6 Washcloths.	50
TT1126	NICU Basic Bath Accessory Pack Each pack contains: 1 Disposable Liner; 1 Scalp brush & 6 Washcloths.	50
TT1105	TurtleTub Disposable Bath Liner - Individually Packaged	100
Swaddle Bathing Blankets		
TT1103	Fleece Reusable Swaddle Bathing Blanket - Small (20" x 23")	25
TT1110	Fleece Reusable Swaddle Bathing Blanket - Large (20" x 30")	25
TT1115	Cotton Flannel Reusable Swaddle Bathing Blanket - Large (20" x 30")	25
TT1119	Non-woven Disposable Swaddle Bathing Blanket - Large (20" x 30")	100
Mobile Bathing Cart		
TT1114	Stainless Steel 2-drawer Bathing Cart - 20" W x 16" D x 34" H, 4-sided guard rail adds 3.5" to Height	1
TT1120	Cart Caddy - 14.5" x 3.75" x 3.75"	1
Baby Wash, Brush & Comb Set, Storage Bag & Rinse Cup		
TT1102	Baby Wash Single Use Packets 10mL Neutral pH Mild Baby Wash	100
TT1127	Brush & Comb Set	100
TT1108	TurtleTub Storage Bags	10
TT1109	TurtleTub Rinse Cups - 22 oz	25

REFERENCES

- 1. Edraki M, et al. Comparing the effects of swaddled and conventional bathing methods on body temperature and crying duration in premature infants: a randomized clinical trial. Journal of Caring Sciences. 2014; 3 (2): 83-91.
- 2. Caka SY, Gözen D. Effects of swaddled and traditional tub bathing methods on crying and physiological responses of newborns. J Spec Pediatr Nurs. 2018;23:e12202
- 3. Swapna, G., et al. Relative Effectiveness of Swaddle Bath and Conventional Bath on Level of Thermal Stability and Crying Duration among Preterm Infants at Selected Hospital in North India. ICCRJNR. Jan-Jun 2017, 2(1): 34-54.

 4. Ceylan SS, Bolisik. Effects of swaddled and sponge bathing methods on signs of stress and pain in premature newborns: Implications for Evidence-Based Practice. World Views on Evidence Based Nursing. 2018; 15: 296-303. doi 10.1111/wvn.12299
- 5. Paran M, Edraki M, Montaseti S, Razavi Nejad M, Comparing the Effects of Swaddle and Conventional Bathing Methods on Behavioral Responses in Preterm Neonates JJN. 2016; NOV, 7(4) DOI:10.22038/jin.2016.7778
- 6. Hall K. Practising developmentally supportive care during infant bathing: reducing stress through swaddle bathing. Infant. 2008; 4(6): 198-201.
- 7. Craig, JW, et al. Recommendations for involving the family in developmental care of the NICU baby. Journal of Perinatology. 2015; 35(suppl 1): S5-S8.
- 8. Quiraishy K, et al. A protocol for swaddled bathing in the neonatal intensive care unit. NAINR. 2013; (1): 48-50.
- 9. Association of Women's Health, Obstetric, and Neonatal Nurses (AWHONN). Neonatal Skin Care: Evidence-based Clinical Practice Guideline. 2018; 41-53.
- 10.Fern D, et al, Swaddled bathing in the newborn intensive care unit. Newborn and Infant Nursing Reviews. 2002; 2(1), 3-4.



Manufactured by:

CATAPULT PRODUCTS, LLC

4677 S. Cherry St. Salt Lake City, UT 84123 385-212-4933 (phone) 385-234-4850 (fax) info@catapultmed.com www.catapult-products.com

